

# UNIVERSAL SURGICAL KIT - VOO DOO TACTICAL UNIVERSAL SURGICAL KIT, COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

You'll always be ready to prevent an emergency from becoming a crisis with our Universal Surgical Kit. Complete with everything you need for quick surgery in the field, with straps to hold all the necessary tools organized and in place, you'll save valuable time when it's needed most. The U.S. Military field style kit contains all new components neatly folded into one convenient kit that attaches to your belt, vest, or gear with a universal strap on its back. Your kit includes: scissors, probe, sutures, scalpel, hemostat, pen light, tweezers, alcohol pads and ointment. the Complete Universal Surgical Kit is constructed of durable G.I. style pack cloth, measures a compact 7"L x 5"W x 2.5"H, and at just over a half-pound, you can be prepared while still traveling light. Inner straps hold essentials in place Universal strap attachment G.I. Style pack cloth Contains: Scissors Probe Sutures Scalpel Hemostat Zipper neck Pen light



## Attributes

- Name: [VOO DOO TACTICAL UNIVERSAL SURGICAL KIT, COYOTE](#)
- Manufacturer: [VOO DOO TACTICAL](#)
- Product no.: 100022581
- Mfr. No.: 10-7688007000
- Size: 1 Person
- Delivery weight: 0.272kg
- Shipping height: 25mm
- Shipping width: 102mm
- Shipping length: 203mm
- UPC: 783377049178

# Table of Contents

- [Startpage](#)
- [Universal Surgical Kit Safety Instructions](#)
- [About Us](#)

# Universal Surgical Kit Safety Instructions

## Introduction

Thank you for choosing the Universal Surgical Kit by VOO DOO TACTICAL. This kit is designed to provide you with essential tools for quick surgery in emergency situations. To ensure safe and effective use of the kit, please read and follow these safety instructions carefully.

## General Safety Guidelines

- Ensure that all components of the kit are intact and in good condition before use.
- Always use the kit in accordance with the instructions provided.
- Keep the kit out of reach of children to avoid accidental injury.
- Be aware of your surroundings and ensure you have a safe environment to perform any procedures.
- If you are unsure about how to use any tool, seek professional guidance.

## Specific Safety Precautions for Use

- **Scissors:** Use caution when handling scissors. Always cut away from your body and keep fingers clear of the cutting path.
- **Scalpel:** The scalpel blade is sharp. Handle it with care and avoid direct contact with skin. Always use the protective cover when not in use.
- **Sutures:** Ensure that you have the necessary training to use sutures properly. Improper use can lead to infection or improper healing.
- **Hemostat:** Use the hemostat carefully to clamp blood vessels. Ensure it is properly sterilized before use.
- **Pen Light:** Use the pen light for visibility in lowlight situations. Do not shine the light directly into anyone's eyes.
- **Tweezers:** Handle tweezers carefully to avoid pinching skin. Ensure they are clean before use.
- **Alcohol Pads:** Use alcohol pads to clean tools or skin before procedures. Dispose of used pads properly.
- **Ointment:** Apply ointment only to clean, dry skin. Avoid contact with eyes and mucous membranes.

## Instructions for Installation and Usage

### 1. Kit Preparation:

- Open the kit and check that all items are present and in good condition.
- Familiarize yourself with the location of each tool within the kit.

### 2. Using the Tools:

- Select the appropriate tool for the procedure you intend to perform.
- Ensure that the working area is clean and welllit.
- Follow any specific instructions for the use of each tool as outlined in this guide.

### 3. After Use:

- Clean all tools thoroughly after use to prevent infection.
- Store the kit in a cool, dry place, ensuring that it is secure and not accessible to children.

## Disposal Instructions

- Dispose of any used or expired components, such as alcohol pads and ointment tubes, in accordance with local regulations.
- Do not dispose of sharp objects like the scalpel or scissors in regular trash. Use a designated sharps

container.

- Ensure that the kit is kept intact and functional until it is no longer needed. If the kit is damaged, follow appropriate disposal methods for each component.

## **Contact Information for Further Support**

For any safety inquiries or concerns regarding the Universal Surgical Kit, please contact the appropriate EU-based safety authority. Ensure you have the product details on hand for reference.

By following these safety instructions, you can ensure a safe and effective experience with your Universal Surgical Kit. Always prioritize safety and seek professional assistance when necessary.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)