RAPID HEMOSTATIC Z-FOLD GAUZE - 1"X5' CELOX RAPID RIBBON HEMOSTATIC GAUZE

Time is critical when treating severe injuries in the field and saving time increases the chance of survival for the casualty in an emergency or hostile situation. To be effective in real use, hemostatic dressings need to work fast. The low profile pack means Celox Rapid Ribbon is easily concealed and the reduced size makes it ideal for treating narrow entry wounds such as bullet, blast fragment and stab injuries. Works with just 60 seconds compression. Stops severe bleeding – fast. Rapid action reduces blood loss. Stops hypothermic bleeding as found in severe traumatic injuries. The current generation of hemostatic agents need at least 3 minutes of compression to work. Celox Rapid Ribbon Gauze removes the delay by speeding up packing time and reducing compression time to just 60 seconds while using the proven hemostatic technology of Celox.

Attributes

- Name: 1"X5' CELOX RAPID RIBBON HEMOSTATIC GAUZE
- Manufacturer: CELOX MEDICAL
- Product no.: 100025883
- Mfr. No.: TS-ZR-015
- Size: 1 Person
- Delivery weight: 0.005kg
- UPC: 755702981842

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RAPID HEMOSTATIC ZFOLD GAUZE Safety Instruction Guide

Introduction

Thank you for choosing the RAPID HEMOSTATIC ZFOLD GAUZE 1"X5' CELOX RAPID RIBBON HEMOSTATIC GAUZE. This product is designed for emergency use to control severe bleeding in critical situations. It is essential to follow the safety instructions outlined in this guide to ensure effective use and safety.

General Safety Guidelines

- Ensure the product is suitable for the intended use.
- Always check the packaging for damage before use.
- Store the product in a cool, dry place away from direct sunlight.
- Keep out of reach of children and vulnerable individuals.
- Familiarize yourself with the product and its application before an emergency arises.
- Report any incidents or unsafe conditions to the relevant authorities.

Specific Safety Precautions for Use

- Use the hemostatic gauze only as directed.
- Do not use if the packaging is damaged or if the product is expired.
- Avoid contact with open wounds unless necessary.
- Ensure that the area of application is clean and dry.
- Do not remove the gauze once applied; it is designed to stay in place until professional medical help arrives.
- If allergic reactions occur, discontinue use and seek medical attention.

Instructions for Installation and Usage

1. Preparation:

- Assess the severity of the bleeding and the injury.
- Ensure you have clean hands or wear gloves if available.

2. Application:

- Remove the gauze from its packaging.
- Place the gauze directly onto the wound, ensuring complete coverage.
- Apply firm pressure to the gauze.
- Maintain pressure for a minimum of 60 seconds to allow the hemostatic agent to work effectively.

3. Postapplication:

- After 60 seconds, check for bleeding. If bleeding continues, apply additional gauze as necessary without removing the first application.
- Once bleeding is controlled, do not remove the gauze; secure it in place with a bandage if available.
- Seek professional medical help as soon as possible.

Disposal Instructions

- Dispose of used gauze and packaging in accordance with local regulations.
- Do not flush the gauze down the toilet.
- Ensure that disposed materials are kept out of reach of children and pets.

Contact Information for Further Support

For any inquiries or safety concerns regarding the RAPID HEMOSTATIC ZFOLD GAUZE, please refer to the product packaging for contact details or consult your local distributor.

Conclusion

The RAPID HEMOSTATIC ZFOLD GAUZE is an essential tool for managing severe bleeding in emergencies. By following these safety instructions, you can ensure effective use and contribute to the safety of yourself and others in critical situations. Always prioritize professional medical assistance following the use of this product.

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