

RAPID HEMOSTATIC Z-FOLD GAUZE - CELOX MEDICAL 3"X5' CELOX RAPID HEMOSTATIC Z-FOLD GAUZE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

Time is critical when treating severe injuries in the field and saving time increases the chance of survival for the casualty in an emergency or hostile situation. To be effective in real use, hemostatic dressings need to work fast. The low profile pack means Celox Rapid Ribbon is easily concealed and the reduced size makes it ideal for treating narrow entry wounds such as bullet, blast fragment and stab injuries. Works with just 60 seconds compression. Stops severe bleeding – fast. Rapid action reduces blood loss. Stops hypothermic bleeding as found in severe traumatic injuries. The current generation of hemostatic agents need at least 3 minutes of compression to work. Celox Rapid Ribbon Gauze removes the delay by speeding up packing time and reducing compression time to just 60 seconds while using the proven hemostatic technology of Celox.



Attributes

- Name: [CELOX MEDICAL 3"X5' CELOX RAPID HEMOSTATIC Z-FOLD GAUZE](#)
- Manufacturer: [CELOX MEDICAL](#)
- Product no.: 100025884
- Mfr. No.: TS-ZR-035
- Size: 1 Person
- Delivery weight: 0.005kg
- Shipping height: 25mm
- Shipping width: 114mm
- Shipping length: 140mm
- UPC: 755702981835

Table of Contents

- [Startpage](#)
- [RAPID HEMOSTATIC ZFOLD GAUZE CELOX MEDICAL Safety Instructions](#)
- [About Us](#)

RAPID HEMOSTATIC ZFOLD GAUZE CELOX

MEDICAL Safety Instructions

Introduction

Thank you for choosing the RAPID HEMOSTATIC ZFOLD GAUZE CELOX MEDICAL. This product is designed to effectively manage severe bleeding in emergency situations. It is crucial to understand the safety guidelines and proper usage instructions to ensure effective treatment and minimize risks.

General Safety Guidelines

- Ensure that the product is stored in a cool, dry place away from direct sunlight.
- Check the packaging for any signs of damage or tampering before use.
- Keep out of reach of children to prevent accidental ingestion or misuse.
- Always read and follow the instructions provided with the product.
- Be aware of potential allergies to materials used in the gauze.

Specific Safety Precautions for Use

- Only use the gauze on external wounds. Do not apply it to internal injuries.
- Ensure the area around the wound is clean before application to reduce the risk of infection.
- Apply firm pressure as directed to achieve effective hemostatic action.
- Do not remove the gauze once applied; it is designed to stay in place until professional medical help is available.
- If bleeding does not stop after application, seek immediate medical assistance.
- Avoid using the product on individuals with known allergies to the components of the gauze.

Instructions for Installation and Usage

1. Preparation:

- Ensure your hands are clean. If possible, wear gloves to prevent infection.
- Open the packaging carefully to avoid damaging the gauze.

2. Application:

- Place the gauze directly over the wound.
- Apply firm pressure for a minimum of 60 seconds to activate the hemostatic properties.
- Ensure that the gauze is in contact with the wound and is not layered over itself.

3. PostApplication:

- After 60 seconds, check to see if the bleeding has stopped.
- If bleeding persists, maintain pressure and seek medical assistance.
- Do not attempt to remove the gauze; allow medical professionals to assess the wound.

4. Storage:

- Store the gauze in its original packaging until needed.
- Ensure it is kept in a designated first aid kit or similar location for easy access.

Disposal Instructions

- Dispose of the used gauze in accordance with local regulations for biohazard waste.
- Do not throw the used gauze in regular trash to prevent exposure to infectious materials.
- Ensure that the disposal method protects others from potential harm.

Contact Information for Further Support

For any inquiries or further assistance regarding the RAPID HEMOSTATIC ZFOLD GAUZE CELOX MEDICAL, please refer to the contact details provided on the packaging or the manufacturer's website.

Conclusion

The RAPID HEMOSTATIC ZFOLD GAUZE CELOX MEDICAL is an essential tool for managing severe bleeding in emergency situations. By following these safety instructions and guidelines, you can ensure effective and safe use of this product. Always prioritize safety and seek professional medical assistance when necessary.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk