

# FOXSEAL CHEST SEAL OCCLUSIVE DRESSING - CELOX MEDICAL

## FOXSEAL CHEST SEAL OCCLUSIVE DRESSING 2/PACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

Foxseal Chest Seal is a pioneering occlusive dressing for the treatment of open chest wounds as a tool in prevention of pneumothorax. Compact, fits IFAK without folding, to avoid creases. Strong adhesion in extreme conditions. Two per pack for ballistic wounds. Simple no-valve design. Vent with large, easy to operate tab. Other dressing options are often folded when carried in kits, leading to a permanent crease which can prevent a true seal. Foxseal is a compact size to fit neatly into pouches, to avoid creasing. Two seals per pack allow entry and exit wounds to be treated. Patient monitoring is vital in these injury types and the design is a simple valve-free option, for simple use. The chest can be vented with the easily-gripped tab and the adhesive re-seals.



### Attributes

- Name: [CELOX MEDICAL FOXSEAL CHEST SEAL OCCLUSIVE DRESSING 2/PACK](#)
- Manufacturer: [CELOX MEDICAL](#)
- Product no.: 100025890
- Mfr. No.: TS-FSOD1
- Size: 2 Person
- Delivery weight: 0.005kg
- Shipping height: 25mm
- Shipping width: 152mm
- Shipping length: 152mm
- UPC: 755702981873

# Table of Contents

- [Startpage](#)
- [FOXSEAL CHEST SEAL OCCLUSIVE DRESSING SAFETY INSTRUCTION GUIDE](#)
- [About Us](#)

# FOXSEAL CHEST SEAL OCCLUSIVE DRESSING

## SAFETY INSTRUCTION GUIDE

### Introduction

Thank you for choosing the FOXSEAL Chest Seal Occlusive Dressing by CELOX MEDICAL. This product is designed to provide effective treatment for open chest wounds and to prevent pneumothorax. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

### General Safety Guidelines

- Ensure the product is stored in a cool, dry place away from direct sunlight.
- Check the packaging for any signs of damage before use. Do not use if the packaging is compromised.
- Always follow the instructions provided with the product for optimal safety and effectiveness.
- Keep out of reach of children and vulnerable individuals unless under adult supervision.
- In case of an emergency, seek immediate medical assistance.

### Specific Safety Precautions for Use

- The product is intended for external use only. Do not apply to internal wounds.
- Ensure hands are clean before handling the dressing to minimize the risk of infection.
- Do not use the dressing on patients with known allergies to adhesive materials.
- Monitor the patient closely for any signs of respiratory distress or complications after application.
- If the dressing becomes loose or is no longer adhering properly, replace it with a new seal.

### Instructions for Installation and Usage

#### 1. Preparation:

- Ensure the area around the wound is clean and dry.
- If possible, have the patient lie down in a comfortable position to minimize movement.

#### 2. Applying the Chest Seal:

- Open the packaging carefully, avoiding any damage to the dressing.
- Remove the protective backing from the adhesive side of the chest seal.
- Center the seal over the wound, ensuring complete coverage.
- Press firmly around the edges to ensure strong adhesion, especially in extreme conditions.

#### 3. Ventilation:

- Use the large, easy-to-operate tab to vent the chest if necessary.
- Ensure the vent is clear and functioning properly to allow air to escape.

#### 4. Monitoring:

- Continuously monitor the patient's condition. Look for any signs of worsening symptoms or complications.
- If the patient shows signs of distress, seek medical assistance immediately.

### Disposal Instructions

- Dispose of used dressing and packaging in accordance with local regulations for medical waste.

- Do not reuse the dressing. Each seal is designed for single use only.

## **Contact Information for Further Support**

For any inquiries regarding safety or product usage, please refer to the manufacturer's contact details provided on the packaging or visit the official website of CELOX MEDICAL.

## **Conclusion**

The FOXSEAL Chest Seal Occlusive Dressing is a vital tool for emergency care in treating open chest wounds. By following these safety instructions, you can ensure effective use and enhance patient safety. Always stay informed about product recalls and safety updates through official channels. Thank you for prioritizing safety and care.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)