RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, ORANGE

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10's products are designed, manufactured and assembled in the USA.

Attributes

- Name: RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, ORANGE
- Manufacturer: ELEVEN 10
- Product no.: 100043728
- Mfr. No.: E10-3011M-ORG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668397

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Safety Instruction Guide for RIGID TQ Case for SOF®TT/SOF®TTW

Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This guide provides important safety instructions and usage guidelines to ensure a safe and effective experience with your product. Please read this guide carefully before using the RIGID TQ Case.

General Safety Guidelines

- Ensure that the RIGID TQ Case is used only for its intended purpose: to securely hold and provide rapid access to the SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of wear, damage, or malfunction. If any issues are detected, discontinue use and replace the case.
- Keep the case out of reach of children and vulnerable individuals who may not understand its purpose or the potential hazards associated with it.
- Always follow local regulations and guidelines concerning the use of medical devices and first aid equipment.

Specific Safety Precautions for Use

- **Avoid Excessive Force:** Do not apply excessive force when inserting or removing the tourniquet from the case to prevent damage to the case or the tourniquet.
- **Wet Conditions:** If the case or tourniquet becomes wet, ensure that they are dried properly before use to maintain effectiveness.
- **Secure Mounting:** Ensure that the case is securely mounted to your duty belt, vest, or holster to prevent accidental detachment during use.
- **OneHanded Deployment:** Practice the onehanded deployment technique to ensure you can access the tourniquet quickly and efficiently in an emergency situation.
- **Friction Mechanism:** Understand that the case utilizes friction to hold the tourniquet in place. Ensure that the case is properly adjusted to hold the tourniquet securely while allowing for quick release.

Instructions for Installation and Usage

1. Mounting the Case:

- Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- Use the appropriate mounting hardware (sold separately) to secure the case in the desired position.
- For Police Duty Belt, ensure that the BladeTech TEKLOK Belt Attachment fits securely onto the belt (up to 2.25 inches).
- For MOLLE/PALS Attachment, ensure that the Tactical Tailor Malice Clip is properly secured.

2. Inserting the Tourniquet:

- Open the case and carefully insert the SOF®TT or SOF®TTW tourniquet into the designated slot.
- Ensure that the tourniquet is seated properly within the case.

3. Using the Tourniquet:

- In an emergency situation, locate the exposed portion of the tourniquet.
- Grasp the exposed portion firmly and pull swiftly to release the tourniquet from the case.
- Follow proper first aid procedures when applying the tourniquet.

4. PostUse Care:

- After use, clean the case and tourniquet as per the manufacturer's instructions.
- Store the case in a dry, safe place away from direct sunlight and extreme temperatures.

Disposal Instructions

- Dispose of the RIGID TQ Case in accordance with local regulations regarding plastic and medical waste.
- Ensure that any used tourniquets are disposed of properly as per local medical waste guidelines.

Contact Information for Further Support

For any inquiries regarding the RIGID TQ Case, please reach out to the manufacturer or your local distributor. It is recommended to keep a record of your purchase for warranty and support purposes.

By following these safety instructions, you can ensure the effective and safe use of the RIGID TQ Case for the SOF®TT/SOF®TTW. Thank you for your attention to safety and for choosing our product.

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