ETAK POUCH MEDICAL KIT CONTENTS - ETAK KIT CONTENTS W/COMPRESSED GAUZE

All the medical items to fill the ETAK Pouch or to replace expired products The Eleven 10 ETAK (Enhanced Trauma Aid Kit) is designed to hold critical medical items with the first responder and tactical officer in mind. Compact and durable, made with premium, duty-proven materials. For use with either the Fixed or Tear Away ETAK Pouches (sold separately). Kit Includes: (1) 4" Compact Medical Trauma Dressing (1) Compression Gauge or Combat Gauze (1) Nitrile Exam Gloves – Pair, Size L (1) Eleven 10® Frog Tape (1) NAR Hyfin Vent Compact Chest Seal, Twin Pack (1) 4" ChemLight Light Stick (1) 5.5" Trauma Shears

Attributes

- Name: ETAK KIT CONTENTS W/COMPRESSED GAUZE
- Manufacturer: ELEVEN 10
- Product no.: 100043774
- Mfr. No.: E10ETAKGA
- Size: 1 Person
- Delivery weight: 0.272kg
- Shipping height: 51mm
- Shipping width: 178mm
- Shipping length: 203mm
- UPC: 680306732968

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ETAK Pouch Medical Kit Safety Instructions

Introduction

Thank you for choosing the ETAK Pouch Medical Kit. This kit is designed to provide essential medical supplies for first responders and tactical officers. Please read this safety instruction guide carefully to ensure the safe and effective use of the products contained within the kit.

General Safety Guidelines

- Ensure safe handling and usage of all medical items.
- Always check the expiration dates of items before use.
- Keep the kit out of reach of children and vulnerable individuals.
- Report any unsafe products or accidents to the relevant authorities.
- Stay informed about product recalls through the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Use only as directed and for the intended purpose.
- Do not use any items that appear damaged or degraded.
- Wear appropriate personal protective equipment (PPE) when using the kit, especially when handling blood or bodily fluids.
- Ensure that you are trained in first aid and the proper use of each item in the kit.
- Do not attempt to use items if you are unsure of their application.

Instructions for Installation and Usage

1. Familiarize Yourself with the Kit Contents:

- (1) 4" Compact Medical Trauma Dressing
- (1) Compression Gauge or Combat Gauze
- (1) Nitrile Exam Gloves – Pair, Size L
- (1) Eleven 10® Frog Tape
- (1) NAR Hyfin Vent Compact Chest Seal, Twin Pack
- (1) 4" ChemLight Light Stick
- (1) 5.5" Trauma Shears

2. Using the Items:

- **Trauma Dressing:** Open the package, and apply directly to the wound, ensuring it covers the area completely.
- **Compression Gauge/Combat Gauze:** Use for severe bleeding. Pack the gauze into the wound and apply pressure.
- **Nitrile Exam Gloves:** Wear before handling any medical items or when treating wounds to protect yourself and the patient.
- **Frog Tape:** Use to secure dressings or other items as necessary.
- **Chest Seal:** Apply to an open chest wound to prevent air from entering the chest cavity.
- **ChemLight Light Stick:** Activate in lowlight situations to signal for help or to illuminate the area.
- **Trauma Shears:** Use to cut clothing away from a wound or to access the injury site.

3. PostUse:

- Dispose of used medical items according to local regulations.
- Replace any expired or used items in the kit promptly to ensure readiness.

Disposal Instructions

- Dispose of used medical supplies in accordance with local waste management regulations.
- Follow specific guidelines for the disposal of biohazard materials if applicable.
- Ensure that the kit is restocked with new items to maintain its effectiveness.

Contact Information for Further Support

For any safety inquiries or further assistance regarding the ETAK Pouch Medical Kit, please refer to the EUbased contact point provided with your purchase.

Thank you for prioritizing safety and preparedness with the ETAK Pouch Medical Kit.

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