

ETAK POUCH MEDICAL KIT CONTENTS - ELEVEN 10 ETAK KIT CONTENTS W/COMPRESSED GAUZE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

All the medical items to fill the ETAK Pouch or to replace expired products The Eleven 10 ETAK (Enhanced Trauma Aid Kit) is designed to hold critical medical items with the first responder and tactical officer in mind. Compact and durable, made with premium, duty-proven materials. For use with either the Fixed or Tear Away ETAK Pouches (sold separately). Kit Includes: (1) 4" Compact Medical Trauma Dressing (1) Compression Gauge or Combat Gauze (1) Nitrile Exam Gloves – Pair, Size L (1) Eleven 10® Frog Tape (1) NAR Hyfin Vent Compact Chest Seal, Twin Pack (1) 4" ChemLight Light Stick (1) 5.5" Trauma Shears



Attributes

- Name: [ELEVEN 10 ETAK KIT CONTENTS W/COMPRESSED GAUZE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043774
- Mfr. No.: E10-ETAK-GA
- Size: 1 Person
- Delivery weight: 0.272kg
- Shipping height: 51mm
- Shipping width: 178mm
- Shipping length: 203mm
- UPC: 680306732968

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [ETAK Pouch Medical Kit Safety Instructions](#)
- [About Us](#)

ETAK Pouch Medical Kit Safety Instructions

Introduction

Thank you for choosing the Eleven 10 ETAK Pouch Medical Kit. This kit is designed to provide essential medical supplies for first responders and tactical officers. To ensure your safety and the effective use of this kit, please read the following safety instructions carefully.

General Safety Guidelines

- Ensure that all items in the kit are used as intended and according to the instructions provided.
- Regularly check the contents of the kit for expiration dates and replace any expired items immediately.
- Store the kit in a cool, dry place away from direct sunlight and moisture.
- Keep this kit out of reach of children and vulnerable individuals to prevent accidental misuse.
- In case of an emergency, always consult a medical professional.

Specific Safety Precautions for Use

- Always wear nitrile exam gloves when handling blood or bodily fluids to reduce the risk of infection.
- Use trauma dressings and gauze only as directed. Improper use can lead to complications.
- Ensure that the compression gauze is applied correctly to control bleeding. Seek medical assistance if bleeding does not stop.
- The ChemLight Light Stick should be used in lowlight situations only. Avoid direct contact with the chemical inside.
- Use trauma shears only for cutting medical tape, clothing, or other materials. Keep away from skin to prevent accidental cuts.

Instructions for Installation and Usage

1. Preparation:

- Open the ETAK Pouch and inspect all contents.
- Familiarize yourself with each item and its intended use.

2. Using the Medical Trauma Dressing:

- Remove the dressing from its packaging.
- Place the dressing over the wound and secure it in place with the provided tape or bandage.

3. Applying Compression Gauze:

- Place the gauze directly over the bleeding area.
- Apply firm pressure and wrap the gauze around the wound to maintain pressure.

4. Wearing Nitrile Exam Gloves:

- Open the glove package and put on the gloves before handling any medical items.
- Dispose of gloves properly after use to prevent contamination.

5. Using the Trauma Shears:

- Hold the shears by the handles and cut away clothing or other materials obstructing access to a wound.
- Ensure the pointed end is facing away from the skin to prevent injuries.

6. Using the ChemLight Light Stick:

- Activate the light stick by bending it until you hear a snap.
- Shake gently to distribute the chemicals and produce light.
- Dispose of the ChemLight properly after use.

Disposal Instructions

- Dispose of expired or used items in accordance with local regulations for hazardous waste.
- Do not throw away items in regular household waste if they are contaminated with bodily fluids.
- Recycle packaging materials where possible.

Contact Information for Further Support

For safety inquiries or further assistance, please refer to the contact information provided with your product packaging.

Thank you for prioritizing safety while using the Eleven 10 ETAK Pouch Medical Kit. Your careful adherence to these guidelines will help ensure effective use and safety in emergency situations.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk