

# TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR

## TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3 x 3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. ESSENTIALS SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Compressed Gauze (1x) SWAT-T Pressure Device (1x)



### Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED BLACK](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050358
- Mfr. No.:
- Size: 1 Person
- Delivery weight: 0.503kg
- Shipping height: 51mm
- Shipping width: 165mm
- Shipping length: 191mm
- UPC: 810073650206

### Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Safety Instruction Guide for Trauma Kit NOW! Essentials Supplies](#)
- [About Us](#)

# Safety Instruction Guide for Trauma Kit NOW! Essentials Supplies

## Introduction

Thank you for choosing the Trauma Kit NOW! Essentials Supplies by Blue Force Gear. This medical kit is designed for quick deployment in emergency situations. It is essential to understand the safety guidelines and proper usage to ensure effective and safe operation. This guide complies with the EU General Product Safety Regulation (GPSR) to help you use this product safely.

## General Safety Guidelines

- Ensure that all users are properly trained in first aid and emergency response procedures.
- Always inspect the kit before use to ensure all components are intact and within their expiration dates.
- Store the kit in a cool, dry place away from direct sunlight to maintain the integrity of the supplies.
- Keep the kit out of reach of children and individuals who may misuse its contents.
- Regularly review the contents of the kit and replace any used or expired items.
- In case of an emergency, always call for professional medical assistance in addition to using the kit.

## Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear gloves when treating injuries to prevent contamination and protect yourself from bloodborne pathogens.
- **Quikclot Bleeding Control Dressing:** Use this dressing only on severe bleeding. Apply direct pressure to the wound and follow the instructions provided with the dressing.
- **Emergency Trauma Dressing:** This dressing is designed for quick application. Ensure it is placed directly over the wound and secured tightly.
- **Hypothermia Blanket:** Use this blanket to maintain body temperature in cases of shock or exposure. Wrap the patient completely to avoid heat loss.
- **Flat Fold Tape:** Ensure the skin is clean and dry before applying tape to avoid irritation.
- **Compressed Gauze:** Use compressed gauze to pack wounds. Always apply direct pressure to control bleeding.
- **SWATT Pressure Device:** Follow the instructions carefully to apply this device properly for effective pressure control.

## Instructions for Installation and Usage

### 1. Mounting the Kit:

- The Trauma Kit NOW! is designed to be mounted on a 3 x 3 MOLLE field or battle belt.
- Use the patented BFG Helium Whisper Attachment System for secure and lightweight mounting.
- Ensure that the kit is firmly attached and does not obstruct movement.

### 2. Using the Kit:

- In an emergency, quickly access the kit and identify the necessary supplies.
- Follow the specific safety precautions outlined above for each component.
- Maintain calm and provide assistance to the injured individual while waiting for professional help.

### 3. PostUse:

- After use, restock any used supplies immediately to ensure the kit is ready for future emergencies.
- Dispose of any contaminated items according to local regulations for hazardous waste.

## **Disposal Instructions**

- Dispose of expired or damaged items in accordance with local regulations.
- Do not dispose of medical supplies in regular household waste.
- Contact local waste management authorities for guidance on the proper disposal of medical waste.

## **Contact Information for Further Support**

For more information regarding the Trauma Kit NOW! Essentials Supplies, please refer to the official Blue Force Gear website or contact your local distributor for assistance. Always ensure that you have access to the latest safety updates and product recalls.

Thank you for your attention to these safety guidelines. Your safety and the safety of others depend on your understanding and adherence to these instructions.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)