TRAUMA KIT NOW! ESSENTIALS SUPPLIES - TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED COYOTE

The Trauma Kit NOW! – MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3 x 3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. ESSENTIALS SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Compressed Gauze (1x) SWAT-T Pressure Device (1x)

Attributes

- Name: TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050359
- Mfr. No.: HW-TKN-MED-ESS-
- Size: 1 Person
- Delivery weight: 0.635kg
- Shipping height: 76mm
- Shipping width: 178mm
- Shipping length: 229mm
- UPC: 810073650213

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Deutsch: TRAUMA KIT NOW! ESSENTIALS SUPPLIES Sicherheitshinweise](#)
- [English: Trauma Kit NOW! Essentials Supplies Safety Instruction Guide](#)
- [Français: Guide de Sécurité du Produit : TRAUMA KIT NOW! ESSENTIALS SUPPLIES](#)
- [Italiano: Guida alle Istruzioni di Sicurezza per il Trauma Kit NOW! Essentials Supplies](#)
- [Český: Návod k použití a bezpečnostní pokyny pro TRAUMA KIT NOW! ESSENTIALS SUPPLIES](#)

TRAUMA KIT NOW! ESSENTIALS SUPPLIES

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Trauma Kit NOW! Essentials Supplies Safety

Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! MEDIUM. This MOLLE mounted medical kit is designed for rapid deployment and holds essential first aid supplies. It is crucial to understand how to use this kit safely and effectively. Please read this guide carefully to ensure safe usage and compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure that all users of the Trauma Kit are trained in basic first aid and emergency response techniques.
- Regularly inspect the kit for any expired or damaged supplies. Replace items as necessary.
- Store the kit in a cool, dry place away from direct sunlight to maintain the integrity of the medical supplies.
- Keep the kit out of reach of children unless under adult supervision.
- In case of a medical emergency, call emergency services immediately after administering first aid.

Specific Safety Precautions for Use

- Always wear Trauma Gloves when handling blood or bodily fluids to minimize the risk of infection.
- Be cautious when using sharp items such as scissors or blades included in the kit.
- Ensure that the Quikclot Bleeding Control Dressing is applied according to the instructions to avoid further injury.
- Use the Hypothermia Blanket to maintain body temperature in cases of shock or exposure.
- Follow the instructions for the SWATT Pressure Device carefully to ensure effective use.
- Do not use any expired or damaged supplies from the kit.

Instructions for Installation and Usage

Installation on MOLLE Compatible Surfaces

1. **Identify the Attachment Point:** Locate a 3 x 3 MOLLE field, battle belt, or other compatible surface.
2. **Align the Kit:** Position the Trauma Kit NOW! so that the attachment straps align with the MOLLE webbing.
3. **Secure the Kit:** Use the patented BFG Helium Whisper Attachment System to secure the kit in place. Pull the straps through the MOLLE webbing and tighten them to ensure a snug fit.
4. **Check Stability:** Ensure the kit is securely attached and does not move when pulled.

Usage Instructions

1. **Accessing Supplies:** Open the kit by unfastening the closure. Ensure you have a clear workspace.
2. **Using Supplies:** Follow the specific instructions for each item included in the kit:
 - **Trauma Gloves:** Put on the gloves before handling any medical supplies.
 - **Quikclot Dressing:** Apply directly to the wound and apply pressure for at least 3 minutes.
 - **Emergency Trauma Dressing:** Place over the wound and secure in place.
 - **Hypothermia Blanket:** Wrap around the individual to retain body heat.
 - **SWATT Pressure Device:** Wrap around the limb and apply pressure as needed.
3. **PostUse:** After administering first aid, dispose of used materials properly and replace any items used from the kit.

Disposal Instructions

- Dispose of expired or damaged supplies according to local regulations.
- Used Trauma Gloves and other contaminated materials should be placed in a sealed plastic bag before disposal to prevent contamination.
- Recycle packaging materials where possible.

Contact Information for Further Support

For any safety inquiries or concerns regarding the Trauma Kit NOW!, please refer to the contact information provided on the packaging or the manufacturer's website. It is important to stay informed about product safety updates and recalls.

Conclusion

The Trauma Kit NOW! MEDIUM is an essential tool for emergency preparedness. By following these safety instructions, you can ensure safe and effective use of the kit. Always stay informed about the latest safety practices and ensure that your kit is well maintained and ready for use. Thank you for prioritizing safety with your Trauma Kit NOW!

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