TRAUMA KIT NOW! ESSENTIALS SUPPLIES - TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED MULTICAM

The Trauma Kit NOW! – MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3 x 3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. ESSENTIALS SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Compressed Gauze (1x) SWAT-T Pressure Device (1x)

Attributes

- Name: TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050360
- Mfr. No.: HW-TKN-MED-ESS-MC
- Size: 1 Person
- Delivery weight: 0.553kg
- Shipping height: 102mm
- Shipping width: 203mm
- Shipping length: 305mm
- UPC: 810073650220

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Trauma Kit NOW! Essentials Safety Instructions

Introduction

Thank you for choosing the Trauma Kit NOW! Essentials. This MOLLE mounted medical kit is designed for quick deployment in emergency situations. It is essential to understand how to use this product safely and effectively to ensure your safety and the safety of others.

General Safety Guidelines

- Always read and understand the instructions before using the Trauma Kit NOW! Essentials.
- Ensure that the kit is kept in a secure location, easily accessible in case of emergencies.
- Regularly check the contents of the kit to ensure that all supplies are intact and within their expiration dates.
- Store the kit in a cool, dry place to prevent damage to the medical supplies.
- Keep the kit out of reach of children to prevent accidental use or ingestion of supplies.

Specific Safety Precautions for Use

- Use the Trauma Gloves provided when handling any blood or bodily fluids to prevent contamination.
- When using the Quikclot Bleeding Control Dressing, apply it directly to the wound and apply pressure until bleeding stops.
- The Emergency Trauma Dressing should be used as directed to ensure proper application and effectiveness.
- The Hypothermia Blanket is intended to retain body heat; ensure it is wrapped securely around the individual in need.
- The SWATT Pressure Device should be applied as instructed to control severe bleeding.
- If you are unsure about how to use any item in the kit, seek professional medical assistance immediately.

Instructions for Installation and Usage

Installation

1. **Mounting the Kit:** The Trauma Kit NOW! Essentials is designed to be mounted on a 3 x 3 MOLLE field, a battle belt, or any other MOLLE compatible surface using the BFG Helium Whisper Attachment System.
 - Ensure the attachment system is securely fastened to prevent accidental detachment during use.
 - Check that the kit is positioned for easy access in emergencies.

Usage

1. **Deploying the Kit:** In an emergency, quickly access the kit and remove it from its mounted position.
2. **Using the Supplies:**
 - **Trauma Gloves:** Put on gloves before handling any medical supplies.
 - **Bleeding Control:** Use the Quikclot Dressing or Emergency Trauma Dressing as necessary.
 - **Hypothermia Blanket:** Wrap the individual securely to maintain body temperature.
 - **SWATT Pressure Device:** Apply to control bleeding where necessary.
3. **PostUse:** After use, check the kit for any used or expired supplies and replace them as soon as possible.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Do not dispose of medical supplies in regular household waste.
- Check for local medical waste disposal programs or facilities that accept such items.

Contact Information for Further Support

For any questions or concerns regarding the Trauma Kit NOW! Essentials, please refer to the manufacturer's website or customer support resources for assistance.

Remember, safety is paramount. Always be prepared and knowledgeable about the use of your Trauma Kit NOW! Essentials to ensure effective response in emergencies.

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