MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MNT COYOTE

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MNT COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050368
- Mfr. No.: HW-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650015

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Micro Trauma Kit NOW! Essentials Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This product is designed to provide essential lifesaving medical supplies for everyday carry. Please read this safety instruction guide carefully to ensure safe and effective use of your Micro TKN.

General Safety Guidelines

- Ensure that the Micro TKN is used only for its intended purpose as a first aid kit.
- Keep the kit out of reach of children and vulnerable individuals to prevent accidental use.
- Regularly check the contents of the kit for expiration dates and replace any expired items.
- Store the kit in a cool, dry place away from direct sunlight.
- Familiarize yourself with the contents and their uses before an emergency arises.

Specific Safety Precautions for Use

- Always wear trauma gloves when providing first aid to avoid contamination.
- Use the Quickclot Bleeding Control Dressing only on bleeding wounds; do not apply to nonbleeding areas.
- The Mini Compression Dressing should be applied firmly but not so tightly that it restricts circulation.
- The Hypothermia Blanket is designed to retain body heat; ensure it is used appropriately in emergency situations.
- Flat Fold Tape should be used to secure dressings; ensure it does not cause skin irritation.

Instructions for Installation and Usage

1. Accessing the Kit:

- The Micro TKN is designed for easy access. Use the BLIP pull tabs to open the pouch with one hand.
- Ensure you are in a safe location before accessing the kit.

2. Using the Contents:

- **Trauma Gloves:** Put on the gloves before treating any wounds.
- **Quickclot Bleeding Control Dressing:** Apply directly to the bleeding area and hold pressure until bleeding stops.
- **Mini Compression Dressing:** Place over the wound and secure with Flat Fold Tape if necessary.
- **Hypothermia Blanket:** Wrap the individual in the blanket to conserve body heat.
- **Flat Fold Tape:** Cut or tear to the desired length and use to secure dressings in place.

Disposal Instructions

- Dispose of expired or used items in accordance with local regulations for medical waste.
- The pouch can be cleaned with a damp cloth; however, do not submerge it in water.
- If the pouch is damaged beyond repair, please recycle it according to local recycling guidelines.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW!, please refer to the EU-based contact point provided with your product packaging.

Conclusion

The Micro Trauma Kit NOW! is a vital tool for emergency situations. By following these safety instructions, you can ensure the effective and safe use of the kit. Remember to stay informed about the contents of your kit and maintain it regularly for optimal performance. Thank you for prioritizing safety.

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