MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT ODG

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT ODG
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100745
- Mfr. No.: HW-TKN-MTKN-PRO-OD
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073650084

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Micro Trauma Kit NOW! Pro Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro. This guide provides essential safety instructions to ensure the safe use of your trauma kit. It is designed for everyday carry and contains vital medical supplies for law enforcement professionals, prepared citizens, and hunters. Please read this guide carefully to understand the safe operation and maintenance of your kit.

General Safety Guidelines

- Ensure that the Micro Trauma Kit NOW! Pro is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure that all items are present and within their expiration dates.
- Keep the kit out of reach of children unless under adult supervision.
- Use the kit only for its intended purpose: to provide first aid in emergency situations.
- Be familiar with the contents and their uses before an emergency occurs.

Specific Safety Precautions for Use

- Trauma Gloves: Always wear gloves when providing first aid to prevent contamination.
- Quickclot Combat Gauze: Use only as directed. Do not apply directly to open wounds without proper training.
- SWATT Pressure Device: Ensure you understand how to use this device before an emergency arises.
- Hypothermia Blanket: Use to maintain body temperature in cases of hypothermia. Ensure it is securely wrapped around the person needing assistance.
- Avoid using any supplies that appear damaged or expired.

Instructions for Installation and Usage

1. **Familiarization:** Before carrying the kit, familiarize yourself with its contents and how to use each item.
2. **Attachment:**
 - The kit can be mounted on a belt or MOLLE system.
 - Ensure the pouch is securely attached to avoid accidental detachment during movement.
3. **Deployment:**
 - In an emergency, use one hand or a single finger to pull the BLIP featured pull tabs for quick access to supplies.
 - Practice deploying the kit to ensure you can do so efficiently in a real situation.
4. **After Use:**
 - After using any supplies, replace them as soon as possible to maintain the kit's readiness for future use.
 - Dispose of used items according to local regulations.

Disposal Instructions

- Dispose of any used or expired medical supplies in accordance with local regulations.
- Do not throw away the kit itself. If the kit is damaged, contact the manufacturer for proper disposal instructions.

Contact Information for Further Support

For any inquiries regarding safety or product support, please refer to the contact information provided with your product packaging.

Conclusion

The Micro Trauma Kit NOW! Pro is a crucial tool for emergency situations. By following these safety instructions, you can ensure that you are prepared to provide effective first aid when it is needed most. Always stay informed about the contents of your kit and practice using it to enhance your readiness.

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