MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT ODG

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT ODG
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100746
- Mfr. No.: BTTKNMTKNPRO
- Size: 1 Person
- Delivery weight: 0.34kg
- Shipping height: 279mm
- Shipping width: 89mm
- Shipping length: 152mm
- UPC: 810073650183

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Micro Trauma Kit NOW! Pro Supplies Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies. This kit is designed for everyday carry, providing essential lifesaving medical supplies for law enforcement professionals, prepared citizens, and hunters. Please read this safety instruction guide carefully to ensure the safe use and handling of this product.

General Safety Guidelines

- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items immediately.
- Keep the kit out of reach of children unless under adult supervision.
- Familiarize yourself with the contents of the kit and their uses before an emergency arises.
- Always follow the instructions provided with individual medical supplies included in the kit.

Specific Safety Precautions for Use

- Trauma Gloves: Use gloves when handling any wounds or bodily fluids to prevent contamination.
- Quickclot Combat Gauze: Apply only as directed. Consult a medical professional if unsure about its use.
- SWATT Pressure Device: Ensure proper application to avoid cutting off circulation.
- Hypothermia Blanket: Use as directed to maintain body temperature in emergency situations.
- Always wash your hands before and after using any medical supplies from the kit.

Instructions for Installation and Usage

1. Mounting the Kit:

- Attach the MOLLE or beltmounted pouch securely to your belt or compatible MOLLE system.
- Ensure that the pouch is easily accessible for quick deployment.

2. Using the Kit:

- To access supplies, pull the BLIP pull tabs with one hand or a single finger from either side of the pouch.
- Remove the necessary items carefully, avoiding any unnecessary force that could damage the supplies.

3. PostUse Care:

- After using any medical supplies, dispose of them properly according to local regulations.
- Restock the kit with new supplies as needed to ensure readiness for future emergencies.

Disposal Instructions

- Dispose of used medical supplies in accordance with local waste management regulations.
- For hazardous materials, follow specific disposal guidelines provided on the packaging.
- Ensure that the kit itself is disposed of responsibly if it is no longer needed.

Contact Information for Further Support

For any inquiries regarding the safety and use of the Micro Trauma Kit NOW! Pro Supplies, please contact your local distributor or refer to the product packaging for further details.

By adhering to these safety guidelines and instructions, you can ensure the effective and safe use of the Micro Trauma Kit NOW! Pro Supplies. Always prioritize safety and preparedness in any emergency situation.

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