TRAUMA KIT NOW! ESSENTIALS SUPPLIES - TRAUMA KIT NOW! ESSENTIAL SUPPLIES OD GREEN

The Trauma Kit NOW! – MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3 x 3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. ESSENTIALS SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Compressed Gauze (1x) SWAT-T Pressure Device (1x)

Attributes

- Name: TRAUMA KIT NOW! ESSENTIAL SUPPLIES OD GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100749
- Mfr. No.: HW-TKN-MED-ESS-OD
- Size: 1 Person
- Delivery weight: 0.492kg
- Shipping height: 57mm
- Shipping width: 133mm
- Shipping length: 210mm
- UPC: 810073650237

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Trauma Kit NOW! Essentials Supplies Safety Instructions

Introduction

Thank you for choosing the Trauma Kit NOW! Essentials Supplies. This mediumsized MOLLE mounted medical kit is designed for rapid deployment in emergency situations. It is important to familiarize yourself with the safety guidelines and usage instructions to ensure effective and safe operation of the kit.

General Safety Guidelines

- **Product Safety:** This product is designed to be safe for all users, including those without technical expertise. Always inspect the kit for any visible damage before use.
- **Enhanced Recalls:** In the unlikely event of a product recall, you will be notified through standardized recall notices. Ensure you stay updated on any such notices.
- **Online Shopping Safety:** If purchased online, ensure that the platform complies with safety requirements equivalent to those of physical stores.
- **Vulnerable Groups:** Extra care should be taken when using this kit around children or individuals with limited mobility.
- **EU Contact Point:** For safety inquiries, please refer to the contact information provided with your product.
- **Rapid Alerts:** Stay informed about potential safety issues by checking the EU's Safety Gate platform for updates on unsafe products.

Specific Safety Precautions for Use

- **Inspection Before Use:** Always check the contents of the kit before use. Ensure all items are intact and within their expiration dates.
- **Proper Handling:** Use gloves when handling medical supplies to minimize contamination and infection risk.
- **Storage Conditions:** Store the kit in a cool, dry place, away from direct sunlight and extreme temperatures to maintain the integrity of the medical supplies.
- **Emergency Situations:** Only use the kit in emergency situations. Familiarize yourself with the contents and their uses before an emergency arises.
- **Age Restrictions:** This kit is intended for use by adults or under adult supervision. Ensure that children do not have unsupervised access to the kit.

Instructions for Installation and Usage

1. Mounting the Kit:

- Locate a suitable MOLLE field, battle belt, or other MOLLE compatible surfaces.
- Use the patented BFG Helium Whisper Attachment System to securely attach the kit.
- Ensure that the kit is mounted firmly to prevent it from detaching during movement.

2. Using the Kit:

- In an emergency, quickly access the kit by unfastening the attachment.
- Familiarize yourself with the contents and their specific uses:
 - Trauma Gloves (1x Pair): Wear these before treating wounds.
 - Quikclot Bleeding Control Dressing (1x): Apply to control severe bleeding.
 - Emergency Trauma Dressing 6" (1x): Use for larger wounds.
 - Hypothermia Blanket Green/Silver (1x): Protect the injured person from hypothermia.
 - Flat Fold Tape (1x): Secure dressings in place.
 - Compressed Gauze (1x): Use to pack wounds and control bleeding.
 - SWATT Pressure Device (1x): Apply for additional pressure to control bleeding.

3. PostUse Care:

- After using the kit, restock any used items as soon as possible to ensure the kit is ready for future use.
- Dispose of any singleuse items according to local regulations.

Disposal Instructions

- Dispose of medical supplies in accordance with local regulations.
- Follow any specific disposal instructions provided with individual items, particularly for hazardous materials.
- Ensure that used items are properly sealed and disposed of to prevent contamination or injury.

Contact Information for Further Support

For further support or inquiries regarding the Trauma Kit NOW! Essentials Supplies, please refer to the contact details provided with your product.

Thank you for taking the time to read this safety instruction guide. Your safety and the safety of others depend on the proper use of this kit. Always prioritize safety and preparedness in emergency situations.

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