MICRO TRAUMA KIT NOW! - PLUS+ - MICRO TRAUMA KIT NOW!-PLUS+-MOLLE-ESSENTIAL-RANGER GREEN

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-PLUS+-MOLLE-ESSENTIAL-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101004
- Mfr. No.: HW-TKN-PLUS-ESS-RG
- Size: 1 Person
- Delivery weight: 0.2kg
- Shipping height: 76mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653047

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Micro Trauma Kit NOW! PLUS+ Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! PLUS+. This kit is designed to provide essential medical supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that all users of the Micro Trauma Kit NOW! PLUS+ are trained in basic first aid and emergency response techniques.
- Keep the kit out of reach of children unless they are under adult supervision.
- Regularly check the contents of the kit to ensure all items are intact and within their expiration dates.
- Store the kit in a cool, dry place to maintain the integrity of the medical supplies.
- In case of any allergic reactions or adverse effects, seek medical assistance immediately.
- Familiarize yourself with the contents of the kit and their intended uses before an emergency arises.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Always wear gloves when treating wounds to prevent infection.
 - Dispose of used gloves properly to avoid contamination.
- QuikClot Products:
 - Follow the instructions provided with each QuikClot product carefully.
 - Do not use QuikClot products on individuals with known allergies to any of the components.
- Emergency Trauma Dressing:
 - Ensure the dressing is applied directly to the wound and secured properly.
 - Do not remove the dressing once applied unless directed by a medical professional.
- Hypothermia Blanket:
 - Use the blanket to maintain body temperature in cases of shock or exposure.
 - Avoid placing the blanket directly against open wounds.
- Flat Fold Tape:
 - Use tape to secure dressings and bandages in place.
 - Ensure the tape does not adhere to the wound itself.
- Pressure Device:
 - Apply the pressure device as instructed to control severe bleeding.
 - Monitor the patient closely for signs of shock or distress.
- Chest Seal:
 - Use the chest seal for penetrating chest injuries to prevent air from entering the chest cavity.
 - Ensure the seal is applied correctly to avoid complications.
- Nasopharyngeal Airway:
 - Use this device only if you are trained in its application.
 - Ensure the airway is lubricated before insertion.
- Decompression Needle:
 - Use this device only if you are trained in its application.
 - Ensure proper technique to avoid injury to the patient.

Instructions for Installation and Usage

1. Familiarization:

- Open the Micro Trauma Kit NOW! PLUS+ and review all components.
- Read the instructions for each item included in the kit.

2. Usage:

- In an emergency, assess the situation and determine the necessary supplies.
- Follow the specific safety precautions outlined for each item as you use them.

3. Storage:

- After use, clean and restock the kit as necessary.
- Ensure all items are returned to their designated places within the kit.

4. Transportation:

- Use the MOLLE med kit mount to securely attach the kit to your gear.
- Ensure the kit is easily accessible in case of an emergency.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not dispose of any medical items in regular trash.
- Contact local waste management authorities for guidance on the proper disposal of hazardous materials.

Contact Information for Further Support

For any inquiries or concerns regarding the Micro Trauma Kit NOW! PLUS+, please consult the manufacturer or your retailer for assistance.

Thank you for your attention to these safety instructions. Proper use of the Micro Trauma Kit NOW! PLUS+ can help ensure the safety and wellbeing of yourself and others in emergency situations.

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