# MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO - COYOTE BROWN

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal.

Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

## Attributes

- Name: MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO - COYOTE BROWN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101007
- Mfr. No.: HW-TKN-PLUS-PRO-CB
- Size: 1 Person
- Delivery weight: 0.24kg
- Shipping height: 89mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653078

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# MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide

## Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+. This comprehensive trauma kit is designed to provide essential medical supplies for emergency situations. To ensure the safe and effective use of this product, please read and follow the safety instructions outlined in this guide.

## General Safety Guidelines

- **Product Safety:** This product is designed for safe use in emergency medical situations. Always ensure that the kit is used according to the instructions provided.
- **Enhanced Recalls:** In the unlikely event of a product recall, standardized notices will be issued. Please stay informed about any updates regarding the safety of this product.
- **Online Shopping:** If purchased online, this product is subject to the same safety standards as those sold in physical stores. Ensure that the retailer complies with safety requirements.
- **Special Consumer Focus:** This kit is intended for adult users. Extra care should be taken when used by or around children or vulnerable individuals.
- **EU Contact Point:** For any safety inquiries, please refer to the appropriate EU contact point listed with the product.
- **Rapid Alerts:** Stay informed about any safety issues related to this product through the EU's Safety Gate system.

## Specific Safety Precautions for Use

- **Do not use if damaged:** Inspect the kit for any signs of damage before use. If any components are damaged or missing, do not use the kit.
- **Keep out of reach of children:** This kit contains items that may pose a choking hazard or cause injury. Store it in a secure location.
- **Follow instructions for each item:** Each component of the kit comes with specific instructions. Follow these instructions carefully to ensure proper use.
- **Use gloves when necessary:** Always wear trauma gloves when providing first aid to prevent contamination and protect yourself and the injured person.
- **Be aware of allergies:** Check for any allergies related to the contents of the kit, particularly with dressings and medications.
- **Seek professional medical help:** This kit is for emergency use only. Always seek professional medical assistance when necessary.

## Instructions for Installation and Usage

### Installation

- **MOLLE Attachment:**
  - To attach the kit to a MOLLEcompatible surface, align the attachment straps with the loops on the surface.
  - Secure the straps by threading them through the loops and pulling tight to ensure the kit is firmly attached.

### Usage

- **Opening the Kit:**
  - To access the contents, unzip the pouch carefully. Ensure you are in a safe location to provide first aid.
- **Using Trauma Gloves:**
  - Roll open the package of trauma gloves and put them on before handling any medical supplies.
- **Applying QuikClot Dressing:**
  - Open the dressing package and apply it directly to the wound, applying pressure as needed to control bleeding.
- **Using Emergency Trauma Dressing:**
  - Remove the dressing from its packaging and place it over the wound. Secure it in place using the attached ties.
- **Using the Hypothermia Blanket:**
  - Unfold the blanket and wrap it around the injured person to help maintain body temperature.
- **Applying Flat Fold Tape:**
  - Use the tape to secure dressings in place or to bind wounds as necessary.
- **Using the Pressure Device:**
  - Follow the specific instructions for the pressure device to manage severe bleeding.
- **Using Full Size Chest Seal:**
  - If applicable, follow the instructions for applying the chest seal correctly to prevent air from entering the chest cavity.

## Disposal Instructions

- **Dispose of used items properly:** After use, dispose of any contaminated materials, such as gloves and dressings, in accordance with local regulations for hazardous waste.
- **Check expiration dates:** Regularly check the contents of the kit for expiration dates. Replace any expired items promptly to ensure the kit is always ready for use.
- **Recycle packaging where possible:** Dispose of packaging materials in accordance with local recycling guidelines.

## Contact Information for Further Support

For further support or inquiries regarding the MICRO TRAUMA KIT NOW! PLUS+, please refer to the contact information provided by the retailer or manufacturer.

Thank you for prioritizing safety and preparedness with the MICRO TRAUMA KIT NOW! PLUS+. Stay safe and be prepared for any emergency situation.

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