MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -WOLF GRAY

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal.

Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101010
- Mfr. No.: HWTKNPLUSPRO
- Size: 1 Person
- Delivery weight: 0.503kg
- Shipping height: 51mm
- Shipping width: 203mm
- Shipping length: 254mm
- UPC: 810073653108

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Safety Instruction Guide for MICRO TRAUMA KIT NOW! PLUS+

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+. This kit is designed for emergency situations, providing essential medical supplies to address trauma effectively. To ensure your safety and the effective use of this product, please read and follow the safety instructions outlined in this guide.

General Safety Guidelines

- Ensure that the Micro Trauma Kit NOW! PLUS+ is used only for its intended purpose.
- Keep the kit out of reach of children to prevent accidental misuse.
- Regularly check the contents of the kit for expiration dates and ensure all items are intact.
- Store the kit in a cool, dry place to maintain the integrity of the medical supplies.
- In case of any adverse reactions or injuries, seek medical attention immediately.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling any medical supplies to prevent contamination.
- Use the QuikClot Bleeding Control Dressing only as directed to effectively manage bleeding.
- Ensure that the Emergency Trauma Dressing is applied correctly to avoid complications.
- The Hypothermia Blanket should be used to maintain body temperature; ensure it is wrapped properly around the patient.
- When using the Nasopharyngeal Airway, ensure it is the correct size for the patient to avoid airway obstruction.
- The Decompression Needle should only be used by trained personnel to prevent injury.

Instructions for Installation and Usage

- **Unpacking the Kit:**
 - Carefully open the kit and check all items against the included inventory list.
- **Using the Medical Supplies:**
 - **Trauma Gloves:**
 - Roll out one pair and wear them before handling any medical supplies.
 - **QuikClot Bleeding Control Dressing:**
 - Apply directly to the wound and apply pressure as needed.
 - **Emergency Trauma Dressing:**
 - Place over the wound and secure it in place.
 - **Hypothermia Blanket:**
 - Unfold and wrap around the patient, ensuring full coverage.
 - **Flat Fold Tape:**
 - Use to secure dressings as needed.
 - **Pressure Device:**
 - Follow the instructions provided to apply pressure correctly.
 - **Chest Seal:**
 - Use as directed to seal chest wounds effectively.
 - **Nasopharyngeal Airway:**
 - Lubricate and gently insert into the nostril, ensuring it is positioned correctly.
 - **Decompression Needle:**
 - Only use if trained; follow the specific instructions for insertion.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Used trauma gloves should be discarded in a sealed bag to prevent contamination.
- Check expiration dates on items and dispose of any expired supplies properly.
- For larger items, consult local waste management guidelines for safe disposal.

Contact Information for Further Support

For any inquiries regarding the MICRO TRAUMA KIT NOW! PLUS+, please reach out to the designated contact point provided with your product. Ensure that you have the product details at hand for a more efficient response.

Thank you for your attention to these safety guidelines. Your safety and the safety of others depend on the proper use of this product.

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