TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES-MULTICAM

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES-MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101018
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 279mm
- Shipping width: 229mm
- Shipping length: 76mm
- UPC: 810073653238

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Trauma Kit NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL by Blue Force Gear. This compact medical kit is designed to provide essential supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- **Product Safety:** This kit is designed to ensure safe use in emergency situations. Always check for any visible damage before use.
- **Enhanced Recalls:** Stay informed about any recalls by checking the EU's Safety Gate platform regularly.
- **Online Shopping:** If purchased online, ensure that the seller complies with EU safety regulations.
- **Special Consumer Focus:** This kit contains items that may not be suitable for children. Keep out of reach of minors.
- **EU Contact Point:** For any safety inquiries, please refer to the manufacturer's contact information provided with your product.
- **Rapid Alerts:** Be aware of updates regarding unsafe products through the Safety Gate system.

Specific Safety Precautions for Use

- **Hazard Identification:**
 - Some items in the kit may contain latex, which can cause allergic reactions.
 - Sharp objects such as scissors or needles may cause injury if not handled properly.
- **Avoiding Hazards:**
 - Always check for allergies before using any medical supplies.
 - Use gloves when handling blood or bodily fluids to minimize the risk of infection.
 - Ensure that all items are used according to their intended purpose as outlined in the instructions.

Instructions for Installation and Usage

- **Opening the Kit:**
 - Pull the top of the pouch to access the contents. Ensure you are in a safe and clean environment when doing this.
- **Using the Supplies:**
 - **Trauma Gloves:** Wear gloves before handling any medical supplies. Dispose of them properly after use.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound and apply pressure to control bleeding.
 - **Emergency Trauma Dressing:** Place over the wound and secure it in place.
 - **Hypothermia Blanket:** Wrap the affected person to maintain body heat.
 - **Flat Fold Tape:** Use to secure dressings in place.
 - **Compressed Gauze:** Pack into the wound to control bleeding.
 - **Chest Seal:** Apply to chest wounds as per instructions provided with the seal.
 - **Nasopharyngeal Airway:** Insert gently to maintain an open airway.
 - **Decompression Needle:** Use only if trained and certified; follow proper protocols.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not throw away items that may still be usable or that could be hazardous in regular waste.
- Check for specific disposal instructions for each item in the kit.

Contact Information for Further Support

For safety inquiries or support, please refer to the contact information provided with your product packaging.

Thank you for taking the time to read this safety instruction guide. Your safety and the safety of those around you is our priority.

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