

HEMOSTATIC GRANULES - CELOX MEDICAL HEMOSTATIC GRANULES, 15G PACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

Celox™ Granules can stop severe bleeding fast. The Celox™ Granules mold to the shape of the wound to get firm pressure on to the source of the bleeding even in complex shaped injuries. Safe, easy to use, cost effective and capable of applying to a variety of different wounds No Exothermic reaction - Does not generate heat Works in hypothermic conditions and clots heparinized blood Simple operation - Tear open, pour, press, pack Celox™ has been proven to work on moderate to severe bleeding and has been extensively tested. Celox granules are available in a 15g sachet for situations where every gram of weight saved is vital, still with enough material to treat a severe injury.



Attributes

- Name: [CELOX MEDICAL HEMOSTATIC GRANULES, 15G PACK](#)
- Manufacturer: [CELOX MEDICAL](#)
- Product no.: 100025879
- Mfr. No.: V12090
- Size: 1 Person
- Delivery weight: 0.005kg
- Shipping height: 25mm
- Shipping width: 102mm
- Shipping length: 140mm
- UPC: 755702981620

Table of Contents

- [Startpage](#)
- [Safety Instruction Guide for Celox™ Medical Hemostatic Granules](#)
- [About Us](#)

Safety Instruction Guide for Celox™ Medical Hemostatic Granules

Introduction

Celox™ Medical Hemostatic Granules are designed to provide rapid control of severe bleeding in emergency situations. This guide offers important safety information to ensure the safe and effective use of the product, in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure the product is used only for its intended purpose: to stop severe bleeding.
- Keep out of reach of children to prevent accidental ingestion.
- Handle with care to avoid contamination of the granules.
- Always check the product packaging for expiration dates before use.
- Report any adverse reactions or unsafe product experiences to local authorities.

Specific Safety Precautions for Use

- Do not apply Celox™ Granules to deep puncture wounds or animal bites without professional medical assistance.
- Avoid contact with eyes; in case of contact, rinse immediately with water and seek medical attention.
- If allergic reactions occur, discontinue use and consult a healthcare professional.
- Use only in situations where professional medical help will follow, as further treatment may be necessary.
- Ensure the wound is cleaned as much as possible before applying the granules to enhance effectiveness.

Instructions for Installation and Usage

1. Preparation:

- Assess the severity of the bleeding and ensure that you are in a safe environment.
- Gather necessary supplies, including Celox™ Granules, gloves, and a clean cloth or bandage.

2. Application:

- Tear open the sachet of Celox™ Granules.
- Pour the granules directly onto the bleeding site, ensuring complete coverage of the wound.
- Apply firm pressure on the wound for at least 23 minutes to promote clotting.
- If necessary, pack the granules into the wound to maintain pressure.

3. PostApplication:

- After the bleeding has stopped, cover the wound with a clean bandage.
- Seek professional medical assistance as soon as possible for further evaluation and treatment.

Disposal Instructions

- Dispose of used sachets and any contaminated materials according to local waste management regulations.
- Do not flush down the toilet or dispose of in a manner that could harm the environment.

Contact Information for Further Support

For any safety inquiries or concerns regarding the use of Celox™ Medical Hemostatic Granules, please contact your local health authority or the manufacturer directly.

Conclusion

By following these safety instructions and guidelines, you can ensure the effective and safe use of Celox™ Medical Hemostatic Granules in emergency situations. Always prioritize safety and seek professional medical help when necessary.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk