

MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MNT MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MNT MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050372
- Mfr. No.: BT-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650121

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Essentials Supplies Safety Instruction Guide](#)
- [About Us](#)

Micro Trauma Kit NOW! Essentials Supplies

Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN) by Blue Force Gear. This guide provides essential safety instructions to ensure the safe use of your trauma kit. The Micro TKN is designed for everyday carry and is equipped with lifesaving medical supplies. Please read this guide carefully to understand how to use the kit safely and effectively.

General Safety Guidelines

- Ensure that the Micro TKN is stored in a cool, dry place when not in use.
- Regularly check the contents of the kit for expiration dates and replace any expired items immediately.
- Keep the kit out of reach of children unless supervised by an adult.
- Always use the kit in accordance with the instructions provided to avoid potential hazards.
- Report any unsafe products or accidents to the relevant authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling any medical supplies to minimize the risk of infection.
- Be aware of potential allergic reactions to materials in the kit, such as latex in gloves or adhesive in dressings.
- If you are unsure how to use any of the supplies, seek guidance from a qualified medical professional.
- Use the kit only for its intended purpose; do not attempt to use it for any other applications.
- Ensure that the kit is easily accessible in case of emergencies.

Instructions for Installation and Usage

1. Installation of the Kit:

- Attach the Micro TKN pouch to your belt or MOLLE system using the provided mounting options.
- Ensure that the pouch is securely fastened and does not move during movement.

2. Using the Kit:

- In an emergency, locate the pull tabs labeled "BLIP" on the pouch.
- Pull the tab with one hand or a single finger from either side to deploy the medical supplies.
- Take out the necessary supplies based on the nature of the injury.
- Follow the instructions on each item for proper application (e.g., how to use the Quickclot dressing or compression dressing).

3. After Use:

- Dispose of any used medical supplies in accordance with local regulations.
- Replenish the kit immediately after use to ensure it is ready for future emergencies.

Disposal Instructions

- Dispose of expired or used items in accordance with local waste management regulations.
- Do not throw medical supplies in regular trash; check for designated medical waste disposal options in your area.
- Ensure that used gloves and other potentially contaminated items are disposed of safely to prevent

health risks.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW! or for safety concerns, please contact the designated EUbased support team.

Thank you for your attention to these safety instructions. Your safety and the safety of those around you is paramount.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk