

IFAS INDIVIDUAL FIRST AID SYSTEM POUCHES - GBRS GROUP IFAS INDIVIDUAL FIRST AID SYSTEM POUCH MULTICAM BLACK

[Shooting Accessories](#) > [Holsters & Belt Gear](#) > [Magazine Belt Pouches](#)

GBRS GROUP IFAS INDIVIDUAL FIRST AID SYSTEM POUCHES The GBRS Group IFAS Individual First Aid System Pouch is a versatile, user-friendly medical pouch designed with the end user in mind. It allows for easy rapid deployment in emergency situations while maintaining a streamlined approach. The squadron laminate construction with tweave stretch pocket for medical kit storage allows for a variety of options. Heavy poly woven elastic retention for the integrated tourniquet holder to colocate your medical gear. FEATURES: MOLLE / BELT compatible Adjustable side release buckle retention 4" Integrated tourniquet sleeve made out of heavy poly woven elastic (fits most types of tourniquets) Laser cut Squadron laminate construction Tweave stretch pocket DIMENSIONS: Main Body - 7.25" x 4.5" Tweave Stretch Pocket - 9.5" x 6" Top Elastic - 4.5" Weight: 2.6oz



Attributes

- Name: [GBRS GROUP IFAS INDIVIDUAL FIRST AID SYSTEM POUCH MULTICAM BLACK](#)
- Manufacturer: [GBRS GROUP](#)
- Product no.: 430114790
- Mfr. No.: GEAR-IFAS-MCB
- Color: -
- Hand: -
- Material: -
- Size: 1 Person
- Delivery weight: 0.68kg
- Shipping height: 51mm
- Shipping width: 184mm
- Shipping length: 305mm
- UPC: 850056753791

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Safety Instruction Guide for IFAS Individual First Aid System Pouches](#)
- [About Us](#)

Safety Instruction Guide for IFAS Individual First Aid System Pouches

Introduction

Thank you for choosing the GBRS Group IFAS Individual First Aid System Pouches. This product is designed to provide you with essential medical supplies in emergency situations. To ensure your safety and the effective use of this product, please read this safety instruction guide carefully.

General Safety Guidelines

- Always inspect the pouch and its contents before use to ensure they are in good condition.
- Keep the pouch out of reach of children unless supervised by an adult.
- Ensure that all users are trained in basic first aid and the proper use of the items contained within the pouch.
- Store the pouch in a cool, dry place away from direct sunlight to maintain the integrity of the materials.
- Regularly check for any recalls or safety alerts related to this product on the EU Safety Gate platform.

Specific Safety Precautions for Use

- Do not use the pouch for any purpose other than its intended use as a first aid kit.
- Ensure that the pouch is securely fastened when attached to a belt or MOLLE system to prevent accidental detachment.
- Use caution when accessing sharp instruments or items within the pouch to avoid injury.
- If the pouch or any of its contents are damaged or compromised, discontinue use immediately and replace them as necessary.
- Be aware of personal allergies to any materials or medications included in the pouch.

Instructions for Installation and Usage

Installation

1. Attachment to MOLLE System or Belt:

- Ensure that the MOLLE straps are threaded through the appropriate loops on your gear.
- If using a belt, slide the pouch onto the belt and secure it using the adjustable side release buckle.
- Confirm that the pouch is firmly attached and will not easily detach during movement.

Usage

1. Accessing Medical Supplies:

- Open the pouch using the top zipper or flap, depending on your model.
- Carefully remove the items you need, taking care to avoid contact with sharp objects.
- Follow any specific instructions included with individual medical supplies.

2. Using the Integrated Tourniquet:

- Locate the integrated tourniquet sleeve.
- Remove the tourniquet and follow the instructions provided with it for proper application.
- Ensure that the tourniquet is used only when necessary and by trained individuals.

3. After Use:

- Repack the pouch with any used items replaced as soon as possible to ensure readiness for future

emergencies.

- Clean any reusable items according to the manufacturer's instructions to maintain hygiene.

Disposal Instructions

- Dispose of any expired or used medical supplies in accordance with local regulations for hazardous waste.
- The pouch itself can typically be disposed of with regular waste unless otherwise specified.
- If the pouch is damaged beyond repair, please recycle it if possible, following local recycling guidelines.

Contact Information for Further Support

For any safety inquiries or further support, please refer to the contact information provided on the product packaging or visit the manufacturer's website.

By following these guidelines, you can ensure the safe and effective use of your GBRS Group IFAS Individual First Aid System Pouch. Thank you for prioritizing safety.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk